



**Description:**

Used in clinical solution extraction and injection, could prevent acupuncture accident effectively.

**Product features:**

Just need to push or slightly knock on the safety cap after injection finished, to reach the purpose of protecting needle tip;

Multiple protection, safety cap will not fall off accidentally;

Effective function, simple design helps reduce the cost of safety syringes greatly.

Luer Slip Syringe + Safety Hypodermic Needles (Specification)

• Needle sizes according to customer's requests

Code	Size	With Needles
S6012325	1ml	23G×1 "
S6032325	3ml	23G×1 "
S6052230	5ml	22G×1 1/4 "
S6102140	10ml	21G×1 1/2 "
S6201840	20ml	18G×1 1/2 "

Luer Lock Syringe + Safety Hypodermic Needles (Specification)

• Needle sizes according to customer's requests

Code	Size	With Needles
S7032325	3ml	23G×1 "
S7052230	5ml	22G×1 1/4 "
S7102140	10ml	21G×1 1/2 "
S7201840	20ml	18G×1 1/2 "

## CAP-TYPE SAFETY SYRINGES


Application Method I



Knock the safety cap slightly After injection

Lock the needle tip, safety protection

Application Method II



Push the safety cap after injection

Lock the needle tip, safety protection