

MEGADYNE

MEGASOFT®

Mega Soft® Patient Return Electrodes



The Standard of Electrosurgery

The Mega Soft is simply set on the operating room table, eliminating the small disposable sticky pad that can damage the patient's skin.

MEGA SOFT RETURN ELECTRODES:

- Easy to use
- Eliminates sticky pads
- No damage to patient skin



Mega Soft Universal can be used for any size patient above 350 gms.



Patient Return Electrode



The Universal and Universal Dual patient return electrodes combine many of the best features of the Mega Soft family into one, easy to use patient return electrode.

The Mega Soft Universal & Universal Dual have been specially designed to provide the highest level of patient safety on patients above 350gms (0.8lb), with no maximum patient weight limit. One pad is all you need for all patient sizes.

It can be used either side up so the return cable can easily be positioned away from the surgical site and towards the generator. The Mega Soft Universal's lighter weight and smaller footprint make it more convenient to position in complex procedures.



*Safety and Convenience
in Electrosurgery*



FEATURES & BENEFITS:

- New, lower profile corner mold
- Lower profile provides increased flexibility to use over positioning pads and frames
- Lighter weight (4.5 lbs) makes pad positioning and transport easy
- Easy to warm through and compatible with many patient warming devices
- Available in both a single and dual cord for use with two generators

*If it isn't Green, it isn't a
MEGADYNE® Reusable
Patient Return Electrode*

MEGA SOFT FEATURES & BENEFITS

- **Zero pad-site burns with Mega Soft**—Safer than REM™ in preventing pad-site burns.
- **Confidence.** We provide a warranty and hold harmless for each pad.
- **Easy set-up.** No prepping or searching for an adequate patient pad site. "As simple as making the bed".
- **Kinder to the patient.** No adhesive-related skin irritations.
- **Environment friendly**—eliminates the waste of disposable sticky pads.
- **Compatible** with most major brands of isolated electrosurgical generators.
- **Radiolucent and Latex Free.**

MEGADYNE® Est. 1985
The Electrosurgical Authority®

Mega Soft Return Electrodes



Mega Soft

The Mega Soft combines the innovative Mega Soft patient return electrode technology with an O.R. table pressure reduction pad. The Mega Soft (0830, 0835, 0840) can be used in a variety of surgical procedures and positions providing a safe, easy to use return electrode.



Mega Soft Dual Cord

The Mega Soft Dual Cord is a patient return electrode with two electrosurgical connectors. This allows one patient return electrode to be used with two separate electrosurgical generators at the same time.



Mega Soft Pediatric

The Mega Soft Pediatric is designed for patients between 0.8 & 50 lbs. With the same pressure relieving Akton polymer and no adhesives, delicate skin is protected.

Ordering Information

Description	Catalog Number	Qty Per Box
-------------	----------------	-------------







Mega Soft® Universal & Universal Dual

Mega Soft Universal Patient Return Electrode (36" L x 20" W, 1/8" thickness)	0845	1
Mega Soft Universal Dual Patient Return Electrode (36" L x 20" W, 1/8" thickness)	0846	1

Mega Soft® & Mega Soft Dual Cord

Mega Soft Patient Return Electrode (46" L x 20" W, 1/2" thickness)	0830	1
Mega Soft Dual Cord Patient Return Electrode (46" L x 20" W, 1/2" thickness)	0835	1

Cables & Accessories

Mega Soft Reusable DetachaCable™ Cord, Standard Connector, 2.4m (8.0 ft.)		M2K-01	1
Mega Soft Reusable DetachaCable™ Cord, Standard Connector, 4.4m (14.4 ft.)		M2K-02	1
Mega Soft Reusable DetachaCable™ Cord, Phone Plug Connector, 2.4m (8.0 ft.)		M2K-03	1
Mega Soft Reusable DetachaCable™ Cord, Phone Plug Connector, 4.4m (14.4 ft.)		M2K-04	1
Mega Soft Reusable DetachaCable™ Cord, Extended Phone Plug Connector, 4.4m (14.4 ft.)		M2K-05	1
Mega Soft Reusable DetachaCable™ Cord, Argon Beam Connector, 2.4m (8.0 ft.)		M2K-06	1
Mega Soft Reusable Compatibility Cable, 2.4m (8.0 ft.)		M2K-08	1
Mega Soft Reusable Compatibility Cable, 4.4m (14.4 ft.)		M2K-09	1
Mega Soft Patch Kit		PKIT001	1

Mega Soft® Pediatric

Mega Soft Pediatric Patient Return Electrode (26" L x 12" W, 1/2" thickness)	0840	1
--	------	---

Mega 2000®

Mega 2000 Reusable Patient Return Electrode (36" L x 20" W)	0800	1
Mega 2000 Sheath	0825	30

Megadyne, the exclusive manufacturer of E-Z Clean non-stick PTFE Coated ElectroSurgical Tips, offers a full line of ElectroSurgical Products including:

- Mega Power® ElectroSurgical Generator
- ElectroSurgical Pencils
- PTFE Coated Laparoscopic Electrodes
- All-In-One™ Hand Control
- All-In-One™ Coated Laparoscopic Electrodes
- Lletz Loop Electrodes
- Mega Soft® Patient Return Electrodes
- MegaTIP® Indicator Shaft® Resposable Laparoscopic Electrodes
- Reusable Indicator Shaft Laparoscopic Electrodes
- Disposable Patient Return Electrodes

To order any of our high-quality E-Z Clean Non-Stick PTFE Coated Cautery Tips or our other fine ElectroSurgical Products, contact Megadyne toll-free at **1-800-747-6110**, or your local Megadyne Distributor.

MEGADYNE®, E-Z CLEAN®, IF IT ISN'T GREEN, IT ISN'T MEGADYNE E-Z CLEAN®, ACE BLADE®, MEGA SOFT®, MEGA POWER®, MEGA FINE®, E-Z PEN®, MEGATIP®, INDICATOR SHAFT® and The ElectroSurgical Authority® are trademarks of MEGADYNE MEDICAL PRODUCTS, INC. and are registered in the U.S. and certain other countries. The color green, as used on MEGADYNE'S E-Z CLEAN® brand non-stick electroSurgical tips and MEGA SOFT® brand patient return electrodes, is itself a registered trademark of Megadyne Medical Products, Inc. ALL-IN-ONE, MEGA VAC, ULTRA VAC, and ATTACHA VAC are trademarks of MEGADYNE MEDICAL PRODUCTS, INC. ACE BLADE coated electrodes are manufactured under one or more of the following U.S. Patents: 6039735, 6066137 and 8500727. INDICATOR SHAFT and MEGA TIP electrodes are manufactured under one or more of the following U.S. Patents: 5531743, 5846237, and 6090107. MEGA SOFT return electrodes are manufactured under one or more of the following U.S. Patents: 6053910, 6214000, 6454764, 6544258, 6582424, 6666859, USD707828 and USD709205. E-Z PEN is manufactured under one or more of the following U.S. Patents: 7124932 and 7255263. Other U.S. and foreign patents issued and pending. AKTON is a trademark of Action Products, Inc.

©Copyright 2015 Megadyne Medical Products, Inc.

3000043-01 rev 002 2015-04

MEGADYNE® Est. 1985

The ElectroSurgical Authority®

11506 South State St.
Draper, Utah 84020
Phone (801) 576-9669
Fax (801) 576-9698

1-800-747-6110
www.megadyne.com

